



75°±10° Back-rest 	35°±10° Knee-rest 	450mm Height Adjustment 350-800mm 	16° Reverse Trendelenburg 	16° Trendelenburg 	Cardiac Chair Position 	Applicable for C-arm Examination 	200mm Bed Extension
weight Weighing System 	CPR Instant Release 	0-90° Angle Indicator 	100+60mm Auto-regression 	30° Position Line 	Built-in LED Torch 		

Oasis Series

Y8y

Electric Bed

SKD-C

Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	350mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Advantages

- ★ Double Auto-regression
- ★ Super-low Bed Position 350mm
- ★ PP Bed Platform, Removable & Washable
- ★ Inserted-type PP Head Foot Board, One Step Release With Safe Lock, Patent Design
- ★ PP Safety Guardrail With Two Steps Unlocking System, Patent Design
- ★ Back Rest 30degree Position, Reduce The Chance Of Lung Infection
- ★ Built-in LED Torch
- ★ Applicable For C-arm Examination

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking System	4pcs
<input checked="" type="checkbox"/> ABS Oasis Side Rail Extra Height	1set
<input checked="" type="checkbox"/> ABS Oasis Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> PP Bed Platform, Removable & Washable	9pcs
<input checked="" type="checkbox"/> Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
<input checked="" type="checkbox"/> Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
<input checked="" type="checkbox"/> Ø100mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> ABS Mattress Stopper	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pc
<input checked="" type="checkbox"/> Manual CPR	

Optional

<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Weighing System
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Cardiac Chair Position
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> 5th and 6th Extra Side Rail For Full Protection
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame
<input type="checkbox"/> Linen Holder



360° View


Oasis Series

Y6y

Electric Bed

SKD-C

75°±10° Back-rest	35°±10° Knee-rest	440mm Height Adjustment 360-800mm	CPR Instant Release	0-90° Angle Indicator	200mm Bed Extension	100+60mm Auto-regression	30° Position Line	Built-in LED Torch
-----------------------------	-----------------------------	------------------------------------------------	-------------------------------	---------------------------------	-------------------------------	------------------------------------	-----------------------------	---------------------------



Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking System	4pcs
<input checked="" type="checkbox"/> ABS Oasis Side Rail Extra Height	1set
<input checked="" type="checkbox"/> ABS Oasis Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> PP Bed Platform, Removable & Washable	9pcs
<input checked="" type="checkbox"/> Plastic Socket Inside PP Guardrail For Phone Or Glasses	2pcs
<input checked="" type="checkbox"/> Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
<input checked="" type="checkbox"/> Ø100mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> ABS Mattress Stopper	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pc

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> 5" and 6" Extra Side Rail For Full Protection
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame

Oasis Series

Y6w

Electric Bed

SKD-C

75°±10° Back-rest	35°±10° Knee-rest	440mm Height Adjustment 365-805mm	CPR Instant Release	200mm Bed Extension	100+60mm Auto-regression	Built-in LED Torch
-----------------------------	-----------------------------	------------------------------------------------	-------------------------------	-------------------------------	------------------------------------	---------------------------



Technical Parameters

External Length	2130mm
External Width	1015mm
Bed Height Min.	365mm
Bed Height Max.	805mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> 5" Covered Castors with Brake	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail with ProErgo-Handle™	1set
<input checked="" type="checkbox"/> ABS Oasis Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> PP Bed Platform, Removable & Washable	9pcs
<input checked="" type="checkbox"/> Ø100mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> ABS Mattress Stopper	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pc

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame



75°±10° Back-rest 	20°±10° Knee-rest 	450mm Height Adjustment 500-950mm 	18° Reverse Trendelenburg 	18° Trendelenburg 	weight Weighing System 	CPR Instant Release 	0-90° Angle Indicator
100mm Auto-regression 	200mm Bed Extension 	Cardiac Chair Position 					

X9x

Electric Bed SKD-C

Technical Parameters

External Length	2040-2240mm
External Width	1020mm
Bed Height Min.	500mm
Bed Height Max.	950mm
Mattress Platform	1820-2020×820mm
Safe Working Load	250kg

Optional

- Weighing System With LCD Screen
- Data Storage System
- Monkey Pole

Technical Configuration

<input checked="" type="checkbox"/> Linak Motor	6pcs
<input checked="" type="checkbox"/> 6" Double Side Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Side Rail with Integrated Control Panel	1set
<input checked="" type="checkbox"/> ABS Bed Ends One Step Removable	1set
<input checked="" type="checkbox"/> Ø75mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> 5" Directional Wheel Pedal	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> X-ray Translucent Mattress Platform	1set
<input checked="" type="checkbox"/> Back-up Battery	1set
<input checked="" type="checkbox"/> Night Light	1set
<input checked="" type="checkbox"/> Bed Push Handle for Emergency Transfer	1set
<input checked="" type="checkbox"/> Manual CPR	
<input checked="" type="checkbox"/> Electric CPR	



360° View

<p>75°±10° Back-rest</p>	<p>35°±10° Knee-rest</p>	<p>290mm Height Adjustment 435-725mm</p>	<p>12° Reverse Trendelenburg</p>	<p>12° Trendelenburg</p>	<p>CPR Instant Release</p>	<p>0-90° Angle Indicator</p>	<p>weight Weighing System</p>
-------------------------------------	-------------------------------------	---------------------------------------------------------	---------------------------------------------	-------------------------------------	---------------------------------------	-----------------------------------------	------------------------------------------

Grace Care

D8d

Electric Bed

SKD-C

Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Grace Side Rail with Integrated Control Panel	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Integrated LCD Control Panel	1set
<input checked="" type="checkbox"/> Ø85mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> Trapeze Socket	2pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> Manual CPR	

Optional

<input type="checkbox"/> Electrical CPR
<input type="checkbox"/> Weighing System with LCD Screen
<input type="checkbox"/> Bed Extension
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Cardiac Position
<input type="checkbox"/> ABS Bed Surface
<input type="checkbox"/> X-ray Bed Surface
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Traction Frame
<input type="checkbox"/> Directional Wheel

Grace Care

D6w

Electric Bed
SKD-C

75°±10° Back-rest	35°±10° Knee-rest	310mm Height Adjustment 435-745mm



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> Hand Controller	1pc
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail with ProErgo-Handle™	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs

Optional

<input type="checkbox"/> Manual CPR	
<input type="checkbox"/> Back-up Battery	
<input type="checkbox"/> X-ray Bed Surface	
<input type="checkbox"/> Monkey Pole	
<input type="checkbox"/> Balkan Traction Frame	
<input type="checkbox"/> Directional Wheel	

Grace Care

D3d

Manual Bed
SK-CD1

75°±10° Back-rest	35°±10° Knee-rest	310mm Height Adjustment 435-745mm



Technical Parameters

External Length	2245mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Grace Side Rail Extra Height	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs

Optional

<input type="checkbox"/> Monkey Pole	
<input type="checkbox"/> Balkan Traction Frame	
<input type="checkbox"/> Directional Wheel	

Grace Care

D3w

Manual Bed SK-CD1

75°±10° Back-rest	35°±10° Knee-rest	310mm Height Adjustment 435-745mm
----------------------	----------------------	-----------------------------------------



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail with ProErgo-Handle™	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs

Optional

<input type="checkbox"/> Monkey Pole	
<input type="checkbox"/> Balkan Traction Frame	
<input type="checkbox"/> Directional Wheel	

Grace Care

D2w

Manual Bed SK-CD1

75°±10° Back-rest	35°±10° Knee-rest
----------------------	----------------------



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height	510mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail with ProErgo-Handle™	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Shoes Holder	1set

Optional

<input type="checkbox"/> Monkey Pole	
<input type="checkbox"/> Balkan Traction Frame	
<input type="checkbox"/> Directional Wheel	



<p>75°±10° Back-rest</p>	<p>35°±10° Knee-rest</p>	<p>290mm Height Adjustment 435-725mm</p>	<p>12° Reverse Trendelenburg</p>	<p>12° Trendelenburg</p>	<p>CPR Instant Release</p>	<p>0-90° Angle Indicator</p>
-------------------------------------	-------------------------------------	---------------------------------------------------------	---------------------------------------------	-------------------------------------	---------------------------------------	-----------------------------------------

B8e

Electric Bed SKD-C

Technical Parameters

External Length	2170mm
External Width	1065mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Side Rail Extra Height	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Integrated Control Panel	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Weighing System
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> X-ray Bed Surface
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Traction Frame
<input type="checkbox"/> Directional Wheel

B6e

Electric Bed SKD-C



Technical Parameters

External Length	2170mm
External Width	1060mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration


<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> Hand Controller	1pc
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs
<input checked="" type="checkbox"/> ABS Side Rail Extra Height	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> X-ray Bed Surface
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Traction Frame
<input type="checkbox"/> Directional Wheel

B4w

Manual Bed SK-CD1



Technical Parameters


External Length	2180mm
External Width	960mm
Bed Height Min.	460mm
Bed Height Max.	710mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	4pcs
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail with ProErgo-Handle™	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs

Optional

<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Traction Frame
<input type="checkbox"/> Directional Wheel



A6k

Electric Bed SK-CD1

Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> ABS Mattress Holder	2pcs

Optional

<input type="checkbox"/> Backup Battery
<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Directional Wheel



360° View

A3k Manual Bed SK-CD1




Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Crank	3pcs	<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> ABS Mattress Holder	2pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set		
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

A2k Manual Bed SK-CD1



Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height	490mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Crank	2pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set	<input checked="" type="checkbox"/> ABS Mattress Holder	2pcs
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		

A1k Manual Bed SK-CD1



Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height	500mm
Mattress Platform	1890×830mm
Safe Working Load	135kg

Technical Configuration


<input checked="" type="checkbox"/> Manual Crank	1pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set	<input checked="" type="checkbox"/> ABS Mattress Holder	2pcs
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		



360° View

H6k

Electric Homecare Bed SKD-C



Technical Parameters


External Length	2015mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Aluminium Alloy Foldable Side Rail	1set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Wood Bed Ends	1set	<input checked="" type="checkbox"/> ABS Mattress Holder	2pcs
<input checked="" type="checkbox"/> Wood Bed Side Panel	1set		

GB8e

Electric Orthopedic Traction Bed SKD-C



Technical Parameters

External Length	2250mm
External Width	1060mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1935×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> ABS Side Rail With Angle Indicator	1set	<input checked="" type="checkbox"/> Aluminium Balkan Frame	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

GT3k

Manual Orthopedic Traction Bed SKG-4C



Technical Parameters

External Length	2270mm
External Width	970mm
Bed Height	500mm
Mattress Platform	1925×830mm
Safe Working Load	250kg


Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pc	<input checked="" type="checkbox"/> Aluminium Balkan Frame	1set
<input checked="" type="checkbox"/> Aluminium alloy ABS Head and Foot Board With Safe Lock	1set		
<input checked="" type="checkbox"/> Aluminium Alloy Foldable Side Rail	1set		
<input checked="" type="checkbox"/> IV Pole	1set		

CX2x

Children Bed

SK-C2



Technical Parameters

External Length	1800mm
External Width	900mm
Bed Height	650mm
External Height	1260mm
Mattress Platform	1760×800mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Full Length Side Rail	1set	<input checked="" type="checkbox"/> Buffer Air Spring	4pcs
<input checked="" type="checkbox"/> Plastic Bed Ends with Cartoon Bear	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CR0q

Children Bed

SK-A2



Technical Parameters

External Length	1600mm
External Width	740mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm


Technical Configuration

<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Full Length Side Rail	1set		
<input checked="" type="checkbox"/> Integrated Bed Ends	1set		
<input checked="" type="checkbox"/> IV Pole	1pc		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CR2q

Children Bed

SK-A2




Technical Parameters

External Length	1600mm
External Width	700mm
Bed Height	500mm
External Height	1000mm
Mattress Platform	1580×660mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Full Length Side Rail	1set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Intergrated Bed Ends	1set	<input checked="" type="checkbox"/> Shoes Holder	1pc
<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs		
<input checked="" type="checkbox"/> IV Pole	1pc		

CD3q Children Bed SK-CD1





Technical Parameters

External Length	2030mm
External Width	1045mm
Bed Height Min.	500mm
Bed Height Max.	770mm
Mattress Platform	1705×900mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> 5" Hardshell Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Aluminium Alloy Full Length Side Rail	1set		
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CT8k Electric Children Bed SKD-C




Technical Parameters

External Length	1700mm
External Width	710mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1485×640mm

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> 5" Double Side Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Aluminium Alloy Side Rail with Extra Protection	1set		
<input checked="" type="checkbox"/> Aluminium Alloy Bed Ends	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CT2k Children Bed SK-C2



Technical Parameters

External Length	1800mm
External Width	730mm
Bed Height	500mm
Mattress Platform	1580×660mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> Aluminium Alloy Side Rail with Extra Protection	1set	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Aluminium Alloy Bed Ends	1set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs	<input checked="" type="checkbox"/> Shoes Holder	1pc