



Aeonmed Anesthesia Machine

CE 0123

aeomed
Reliable Quality Thoughtful Service

Aeon8700A

Anesthesia Machine

Multi- Ventilation Modes
(IPPV, PCV, SIMV-VC, SIMV-PC,
PCV-VG, PS, MANUAL, etc.)

Built-in auxiliary
O₂ flowmeter

APL

Integrated CO₂ circle
absorber with heating system

Bright adjustable
endoscope light system

12.1" TFT LCD
with touch screen

Spotlight (Optional)

Auxiliary table (Optional)



- The comprehensive ventilation modes including PCV-VG, SIMV-PC provide more flexibility to ventilation strategy, suitable for various patients operation including pediatric and adult
- Advanced ventilation management, provides accurate delivery and broad of settings enabling effective care across a wide range of patient sizes and acuity types
- Support both low-flow anesthesia and mini-flow anesthesia
- Compact breathing system, easy to remove and clean, fully autoclavable, optional CO₂ bypass
- Electronic flow-meter, precise control and convenient operation
- 12.1" TFT LCD with color touch screen, highlighted display
- Optional comprehensive gas monitoring include: oxygen (paramagnetic), carbon dioxide, AG, etc.
- CE marked, meets EU clinical requirement

Aeon8600A

Anesthesia Machine

- PCV and SIMV+PS modes ensure flexibility in ventilation, machine doable for both pediatric and adult operation
- Optional functions monitoring O_2 , CO_2 and anesthetic gas
- Flowmeter provides minimum 25% concentration of oxygen in O_2/N_2O mixture
- Absorb system sterilized with high temperature and high pressure
- Enhances usability with auxiliary oxygen, suction, auxiliary power source, monitor arm, infusion pump etc.
- CE marked, meets EU clinical requirement



Aeon8300A

Anesthesia Machine

- PCV and SIMV+PS modes ensure flexibility in ventilation, suitable for both pediatric and adult operation
- Optional functions monitoring, CO_2 and anesthetic gas
- Flowmeter provides minimum 25% concentration of oxygen in O_2 / N_2O mixture
- Absorb system easy to remove and clean, fully autoclavable, optional CO_2 bypass
- Enhances usability with auxiliary oxygen, suction, auxiliary power source, arm for monitor, infusion pump etc.
- CE marked, meets EU clinical requirement



Aeon7200A

Full Configuration

- CE marked, meets EU clinical safety requirement
- Large LCD screen, various monitoring waveform and parameters in different colors
- IPPV, Pressure mode, SIMV modes can easily meet the ventilation requirement
- Integrated CO₂ absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast-release APL valve
- 6-tube of flow meters for 3-gas sources meet long operation requirement
- Optional AGSS provides more protection
- Optional multi-gas monitoring
- Optional CO₂ bypass function
- Spirometry loops



Aeon7200A

Basic Configuration

- CE marked, meets EU clinical safety requirement
- Large LCD screen, various monitoring waveform and parameters in different colors
- IPPV, Pressure mode, SIMV modes can easily meet the ventilation requirement
- Integrated CO₂ absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast-release APL valve
- 4-tube of flow meters for 2-gas sources meet long operation requirement
- Optional AGSS provides more protection
- Optional multi-gas monitoring



Aeon7200

Full Configuration

- CE marked, meets EU clinical safety requirement
- Large LCD screen, various monitoring waveform and parameters in different colors
- Automatic and manual modes can easily meet the ventilation requirement
- Integrated CO₂ absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast-release APL valve
- 4-tube of flow meters for 2-gas sources meet long operation requirement
- Optional AGSS provides more protection
- Optional FiO₂ monitoring

**Aeon7200**

Basic Configuration

- CE marked, meets EU clinical safety requirement
- Large LCD screen, various monitoring waveform and parameters in different colors
- Automatic and manual modes can easily meet the ventilation requirement
- Integrated CO₂ absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast-release APL valve
- 2-tube of flow meters for 1-gas source meet long operation requirement
- Optional AGSS provides more protection
- Optional FiO₂ monitoring



VP500

Vaporizer

- Exemplary handling qualities in operation or transport
- New ergonomic design, operation more comfortable
- T-position, any position and angle of tilt
- Enflurane, Isoflurane, Sevoflurane and Halothane available
- Connection type: Selectatec compatible or Draeger plug-in compatible
- Filling type: Pour Fill, Key Fill, Quik-Fil



VP300

Vaporizer

- With CE certificate
- Meet ISO 8835-4 International Standard
- Interlock system
- Pressure, temperature and flow compensation
- Enflurane, Isoflurane, Sevoflurane and Halothane available
- Connection type: Selectatec compatible or Draeger plug-in compatible
- Filling type: Pour Fill, Key Fill, Quik-Fil



AP1000 AGSS

Anesthesia Gas Scavenging System

- With CE certificate
- Meet ISO 8835-3 International Standard
- Great compatibility
- Suction pressure can be adjusted to different requirement
- Good readability with buoy indication
- Satisfy the purification requirement of laminar flow operation room
- Simple operation, easy connection and convenient maintenance



aeomed
BEIJING AEONMED CO., LTD.

Office (Headquarters):
ADD: 11B2, Fengtai Science Park, (100070) Beijing, China
TEL: +86-10-8368 1616 FAX: +86-10-6371 8989
Http://www.aeonmed.com E-mail: marketing@aeonmed.com

Information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy Beijing Aeonmed Co., Ltd. reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.

AAM-1506