

D001 Orthopaedics traction frame (table matching accessories)



D00B1 can be matched up to various type operating tables, accompanied with kinds of orthopedic operating accessory parts, for crotch and lower limbs operations. Its composed of connected fixed device, perinaeum colum, external rod, retractor orthoterion, foot pedal fixed plate, side slide way, support column, steering joint, carframe, traction foot pedal, traction foot pedal universal chuck.

Material use high standerd stainless steel, with smooth design and connection for long time service.

Traction movement	$\geq 900\text{mm}$
Traction frame thigh board adjusted	$\geq 320\text{mm}$
length	$\geq 1650\text{mm}$
height	700 ~ 1020mm
Traction shoe board adjusted	$\geq 350\text{mm}$