


D0010 Carbon fiber upper table (table matching accessories)



Clinically adapt to upper limb (finger, wrist, arms and shoulder) surgeries assistant. Working board adopts full carbon fiber material, can make high clear X-ray imaging, wide range non-metal area have a maximum X-ray imaging effect. Easy disassemble, with a low weight and stable structure. Unique design side cantilever bracket steel board, can adjust height lifting up and down with operating table, to smoothly support the surgery.

Traction movement	≥ 900mm
Traction frame thigh board adjusted length	≥ 320mm
height	≥ 1650mm
Traction shoe board adjusted	700 ~ 1020mm
	≥ 350mm