

DST-II Ophthalmic electric operating table



- 1.this product at present, is one of the most advanced ophthalmologic operating table in domestic.
The lowes height is 600mm, the head board has the eletric adjustment function.
- 2.the head board and table lifting are all controlled by a security-voltage motor and head of the table is arranged with a ,multi-purpose hand braket .

Total length(head board down the lowest)	2000 ± 20mm
Total length(head board up the heighest)	2070± 20mm
width	600 ± 20mm
Height (with mattress)	640 ± 20mm
Head plate lifting route (electric)	60 ± 10mm
Lifting route(electric)	200± 20mm