

DST-IV Obstetrics and gynecology electric delivery table



size	1080×2058×550 mm
Minimum height	550±10mm
Lifting route(electric)	360±10mm
Voltage	AC 230V 2.5A 50Hz
Net weight/ capacity weight	305kg/ 380kg
Side arm-rest support	180° rotation
Rotated tray support board longitude(manual)	60±5mm
leg support board (electric)	left: 0-33 °
	right: 0-33 °
seat board (electric)	Left: ≥25 °
	Right: ≥ 25 °
leg board fold	0-90°(air spring)
back board (electric)	Upward: ≥ 25 °
	Downward: ≥ 25 °
tray support board (manual)	lateral movement: 150+/-3mm
	rotation: 90°+/-3°
function	CPR one reset function
anti-water	IPX0, IPX1
cover	PVC leather

1. This table is used for puerperal labor, gynecological surgery, diagnosis and gynecological examination as well as other medical functions including emergency cesarean section.

2. be equipped with a concealed movable auxiliary table-board, facilitating operations of the gynecologists and midwives.

3. main body adopts 45# steel, cover material adopts cold rolled steel and 304 stainless steel, base cover is ABS with anti-dazzling.

4. Imported linear motor has low noise, smooth and precise. Unique three structure, good heat retaining and cushioning.

5. The lowest height is 550mm for many gynecology, the table board is composed of back, seat and foot board, by handling shank to control back board up and down and rotation, trendlenburg and anti-trendlenburg. Neonatal delivery tray can be rotated for side nurse.