

ATLAS N3

Configurations

Dimensions	Height (With casters)	146cm±1cm
	Width	90cm±1cm
	Depth	71cm±1cm
	Safe shelf loading	25kg±0.5kg
Screen	8 inch TFT LCD screen, 800x480 pixels	
Gas supply	O ₂ , N ₂ O, AIR	
Selectatec bar	Dual positions	
ACGO	Standard	
Mechanical Flowmeter	6 flow tubes for 3 gas	
By-pass	Standard	
Ventilator Software	V-CMV, P-CMV, PSV, V-SIMV, P-SIMV, PRVC, Standby, Manual/Spont	
Spirometry loop	P-V,P-F,F-V,reference Loop	
Spare Cylinder yoke	Optional(O ₂ , N ₂ O)	
Li-ion Battery	1 battery,4800mAh, Optional(2 batteries,9600mAh)	
AGSS	Optional	
Waveforms	Up to 3 waveforms, user configurable	
Auxiliary power outlets	3	
Wheels	Four wheels(125mm) with 2 separate brakes	
Drawers	2	
Reading lamp	LED lighting Included	
Gas monitor module	Optional (CO ₂ , AG)	
Built-in Heater	Standard	
O ₂ cell	Optional	
Vaporizer	Optional (Drager/Penlon/Northern)	
Patient monitor	Optional	
Suction Device	Optional	

For more information, please contact Northern sales representatives

Northern Meditec Limited

ADD: 4th Floor, Building C, Jin Wei Yuan Industrial Park, Julongshan Area, Kengzi St.
Pingshan District, Shenzhen City, P.R.China

TEL: +86 755 29970566 FAX: +86 755 86528647
EMAIL: info@northernmeditec.com WEBSITE: www.northernmeditec.com

