



Northern Meditec

Northern High Flow Nasal Cannula

Crius HF3 is a high-flow nasal cannula oxygen therapy device capable of delivering up to 100% humidified and heated oxygen at a flow rate of up to 80 L/min, it delivers an oxygen mixture with temperature and easily controlled humidity. The supply system Oxygen provides a great solution for airway control and improved oxygenation.

easy and safe

- Touch the user interface, the navigation wheel and the Quick button to provide an easy operating experience.
- Friendly interaction, fast startup.
- Complete alarm and warning system to ensure the safety of patients during therapy.



Crius HF3



suitable departments

Department of Internal Medicine, Surgical Department, ICU, Emergency Department, Department of Pediatrics, Obstetrics and Gynecology, Department of Geriatrics, Department of Rehabilitation, etc.

Indications

- Type 1 respiratory failure .
- Mild ARDS.
- Mild ventilation dysfunction.

Intolerance to traditional oxygen therapy or ventilation non-invasive CPAP.

- Patients need airway protection or humidified airway.

Screen	7 inch touch
non-invasive mode	CPAP, S, S/T, T
FiO2	21%-100%
Temperature	31y-37y HF 34y LF
gas flow	2-80L/min
active humidifier	1-6 Level
SpO2	Optional
Monitor	FiO2, gas flow, temperature, SpO2, pulse, alarm
shortcuts	On-off, silent
trend chart	1, 3, 7 day data review
Snorkel Alarm	Automatic detection of breathing tube Pressure, therapy time, oxygen temperature, etc.

precise treatment

- Accurate and stable O2 concentration output.
- High flow stable gas.
- Stable PEEP for oxygenation and ventilation.
- Reduces the work of breathing.
- Heating and humidification of mixed gases.

Advanced Monitoring

- High sensitivity sensor.
- SpO2 and PI monitoring in real time.
- O2 flow and concentration display in real time.

comfortable treatment

- It is not necessary to establish an artificial airway.
- Without a mask to avoid claustrophobia.
- There is no man-machine confrontation.
- Patients could eat and speak normally to ensure better rehabilitation.



Northern Meditec Limited

Add: 4th Floor, Building C, Jinweiyuan Industrial Park, Julongshan Area, Longtian Street, Pingshan District, Shenzhen, P.R.China

Phone: +86 755 29970566 Fax: +86 755 86528647

Email: info@northernmeditec.com

Website: <http://en.northernmeditec.com>

