

The MEDITOM DT-200S/150S is a multi-purpose electrosurgical generator.
Generator for use in modern operating rooms and surgical centers.

Monopolar and functions bipolar to meet all your .are available surgical demands - 2
cutting modes/COAG monopolar
• 1 bipolar cutting/coagulation mode

The DRM (Daiwha Patient Return Electrode Monitoring) system is designed to disable the generator
to prevent injury or burns during operation.
An alarm will sound and the output will cut off the be automatically invent of
(ONLY APPLICABLE WITH DOUBLE)CUSHION



1. Trastorno bipolar
2. Monopolar

Capable of use with accessories of various brands.
Different user memory slots can be configured for convenience.

Varios modos de funcionamiento

General Surgery
Dermatology
Endoscopy
Vascular Surgery
Gynecology
Cardiac and Thoracic Surgery
Otorhinolaryngology
Pediatric Surgery
Minimally invasive surgery (MIS)
Neurosurgery
Plastic surgery
Orthopedics
Urology Surgery
Surgery
Transurethral resection (TUR)

Product name	Electrosurgical unit
Nombre del modelo	MEDITOM DT-200S/150S
Main frequency	380 380 380 KHz KHz KHz *10% , , , .475 475 KHz KHz *10%*10%*10%*10%
Repetition rate	33 kHz
HFLC	Less than 150 mA
LFLC	-Patient grounding: Less than 0.01 mA -Chassis to ground: Less than 0.1 mA
Power supply	120/230 V (*10 %). 50/60 Hz
	137 x 300 x 395 millimeters
Weight	13 kilos
Power consumption	800 VA

Tolerance: + 20%>

<tolerance:+ 20%>

MODO		PRODUCCION POTENCIA MÁXIMA	Señor calificado (ohmios)	MODO		PRODUCCION POTENCIA MÁXIMA	Señor calificado (ohmios)
MONOPOLAR	PURO	200 W	500'D	MONOPOLAR	PURO	150 W	500 n
	MIX	150 W	300H		MIX	120 W	300 Ct
MONOPOLAR COAG	CONTACT	80 W	3000	MONOPOLAR COAG	CONTACT	70 W	300 n
		70 W	5005			60 W	500 Ct
BIPOLAR		80 W	200 0	BIPOLAR		70 W	200 Q
	COAG	70 W	100 ID		COAG	60 W	100G