

INFUSION PUMP SYRINGE

Ref: M300



General Description:

The Comen M300 Syringe Infusion Pump is designed to deliver high-precision injections with advanced multi-CPU monitoring. This device combines cutting-edge technology, a dynamic pressure monitoring system, and a customizable medication library to ensure maximum safety and comfort in medical procedures. Its high-performance motor and anti-bolus and dynamic monitoring systems make it a reliable choice for a wide range of clinical applications.

Standard Configuration:

- **Main Unit:** M300 infusion pump with advanced CPU monitoring.

Included Components:

- Power cable.
- User manual in local language and English.
- Universal syringe holder.
- **Compatibility:** Compatible with 5ml, 10ml, 20ml, 30ml and 50ml syringes.

Alarms and Security:

- Antibolus system with automatic pressure release.
- Real-time dynamic pressure monitoring.

Visual and audible alarms for:

- Obstructions.
- Low battery.
- Completion of infusion.
- Abnormal pressure.

Integrated Features:

- Customizable library of 2100 medications.
- Continuous titration for uninterrupted adjustments.
- Cascade infusion for administration of multiple drugs or the same drug in sequence.

Power Specifications: Voltage:

- 100-240V AC.
- **Frequency:** 50/60 Hz.
- Long-lasting rechargeable battery for continuous use in case of power outages.