

BEIJING

Infant Incubator Model B1000



STANDARD COMPOSITION:

- Servo Control Mode: air temperature.
- Instantaneous parameter: set temperature, air temperature and heat feed rate.
- Alarms: Audible and visual power outage, sensor failure, major and minor temporary, air temperature override, fan failure, system error.
- >37°C Default mode: power failure total clear, alarm silence, keyboard lock, self-test function.
- Double slope hood: 6 hand ports; left and right cabinets, top platform.
- RS-232 Output connector (use for network system).