

## BEIJING

### Infant Incubator Model B2000



#### STANDARD COMPOSITION:

- Servo Control Mode: air temperature, skin temperature;
- Instantaneous parameter: set temperature, air temperature, skin temperature, heater feed rate;
- Alarms: Audible and visual power outage, sensor failure, major and minor temporary, air override and skin temperature, fan failure, system failure;
- >37°C Default mode: power failure total clear, alarm silence, keyboard lock, self-test function;
- Double slope hood: 6 hand ports; left and right cabinets, top platform;
- RS-232 Output connector (use for network system).