





Infant Incubator BB-200 Luxurious

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Air temperature, skin temperature, time, heating power separately displayed by LED
- Five alarm functions by audible and visual
- $>37^{\circ}\text{C}$ temperature setting
- The infant bed inclination is adjustable
- With phototherapy unit and tray
- Humidity is adjustable continuously
- Optional function: Phototherapy unit down side BL-100; Height adjustable; RS-232 connector

Infant Incubator BB-200 Standard

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Air temperature, skin temperature, time, heating power separately displayed by LED
- Five alarm functions by audible and visual
- $>37^{\circ}\text{C}$ temperature setting
- The infant bed inclination is adjustable
- Humidity is adjustable continuously
- Optional function: Phototherapy unit up side BL-30; Down side BL-100; Access Door; Tray and pole; RS-232 connector; Oxygen inlet



Transport Incubator BT-100

- Air mode and Baby mode controlled by micro-computer
- Use AC Power and DC Power alternatively, DC 12V or DC 24V Power on the ambulance
- Set temperature, air temperature, baby temperature and internal battery voltage can be displayed by LED separately
- Double wall hood with side door the bassinet can be pulled out easily
- Natural air flow humidity
- The height of whole unit can be adjusted
- Oxygen cylinder and Oxygen supply system
- Observation lamp, adjustable brightness

